


# Workflow



# Commands

- Download from SRA (SRATools)

Install SRATools and execute

« *fastq-dump SRRXXXXXXX* » to download SRRXXXXXXX data set (single end), add  
« *--split-files* » for paired-end

- Quality check (FASTQC)

Install FASTQC and execute :

« *fastqc file\_name* » with filename fasta or fastq file

- Trimming (trimmomatic) // cluster

install trimmomatic and execute following command to cut the file at X pb :

« *fastx\_trimmer -l X -i inputfile.fastq -o outputfile.fastq* »

# COMMANDS WITH GTF

- Index creation (STAR)

Execute on cluster :

```
« STAR --runMode genomeGenerate --runThreadN 40 --genomeDir INDEX_DIR --genomeFastaFiles  
REFERENCE_GENOME_FASTA_FILE --sjdbGTFfile ANNOTATION_GENOME_GTF_FILE --sjdbOverhang X »
```

with X the (length-1) of your reads

- mapping (STAR) // cluster

Execute on cluster :

```
« --runThreadN 40 --twopassMode Basic --outSAMstrandField intronMotif --outFilterMultimapNmax -1 --genomeDir  
INDEX_DIR --outSAMattributes All --outSAMtype BAM SortedByCoordinate --outFileNamePrefix OUTPUT_FILE_PREFIX  
--readFilesIn INPUT_FILES_FASTQ »
```

outFilterMultimapNmax : determine the number of multihit (-1 for complete multiple alignment)

- Counts (featureCounts) // cluster

Install subread package and execute :

```
« featureCounts -T 20 -g gene_id -a ANNOTATION_FILE.gtf -o OUTPUT.txt INPUT.sam(.bam) »
```

To count reads per gene with 20 threads from a BAM or SAM file

# COMMANDS WITHOUT GTF (1/2)

- Index creation (STAR) without annotation file // cluster

Execute on cluster :

```
« STAR --runMode genomeGenerate --runThreadN 40 --genomeDir INDEX_DIR --genomeFastaFiles  
REFERENCE_GENOME_FASTA_FILE --sjdbOverhang X »
```

with X the (length-1) of your reads

- Mapping (STAR) without annotation file // cluster

Execute on cluster :

```
« --runThreadN 40 --twopassMode Basic --outFilterMultimapNmax -1 --genomeDir INDEX_DIR  
--outSAMattributes All --outSAMtype BAM SortedByCoordinate --outFileNamePrefix  
OUTPUT_FILE_PREFIX --readFilesIn INPUT_FILES_FASTQ »
```

outFilterMultimapNmax : determine the number of multihit (-1 for complete multiple alignment)

- Bam files sorted By gene/contig (samtools)

```
« samtools sort sample.bam »
```

- Sam files (samtools)

```
« samtools view -h -o out.sam in.bam »
```

# COMMANDS WITHOUT GTF (2/2)

- Configure header and body of the sam file

Manual configuration (with python script or UNIX commands) for custom parseur

- Parsing sam files with custom parser

« *python3 parseur.py* »

The output is 2 files per file.sam, unique alignment count and multiple alignment count

- Fusion of the differents output

With custom python script

« *python3 fusion\_file.py* »

The output is a count table delimited by tab

# COMMANDS

- Organise Count table (R)
- DESeq2 analyses (R, rstudio)
  - « *dds=DESeqDataSetFromMatrix(countData=countdata,colData = coldata,design=~condition)* »
  - « *deseq=DESeq(dds)* »
  - « *res=results(deseq,name="condition\_root\_vs\_nodule")* »
  - « *resSig=subset(res,pvalue<0.01)* »
  - « *write.csv(resSig,file="resSig.txt")* »